---

**CheckMate® Metatarsophalangeal Arthrodesis System**  
**Bone Screws & Surgical Instruments**

**DESCRIPTION**

The CheckMate® Metatarsophalangeal (MTP) Instrument Tray consists of implantable bone screws contained in a screw caddy, and surgical instruments, which are designed to be used in sizing, location, and delivery of the CheckMate® MTP bone plate and bone screw fixation components. The surgical instruments are designed for repeated use, with proper care and handling.

**INDICATIONS FOR USE**

The CheckMate® Metatarsophalangeal (MTP) Arthrodesis System is intended for use in stabilization and fixation of the 1st MTP joint in the foot for fusion, osteotomy, nonunion, malunion or revision surgery.

Instructions for implant delivery using the CheckMate® instruments is provided in the CheckMate® Implant *Instructions for Use (IFU)* package insert, as well as the available *Surgical Technique Guide*.

**WARNINGS AND PRECAUTIONS**

The CheckMate® MTP Instrument Tray is provided **NON-STERILE**. It must be properly cleaned and sterilized before each use. Remove and discard any plastic caps or tip protectors before cleaning and sterilizing the instruments. Bone screws are intended for Single Use only. Any unused bone screw i.e. not implanted, but bought in contact with a patient needs to be discarded and disposed. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing the CheckMate® implantation procedure. The CheckMate® instrument tray and instrument set should be regularly inspected for any signs of wear or damage.

The CheckMate® bone screws (locking, non-locking and interfragmentary) are designed to be used **ONLY** with the CheckMate® bone plates. Use of these screws as standalone or with other bone plates and/or plating systems will compromise patient safety and device performance.

**SCREW CADDY ASSEMBLY**

- The screw caddy forms a part of the instrument tray, and needs to be populated with bone screws adequately and appropriately, before each cleaning and sterilization cycle.

- The bone screws (locking, non-locking and interfragmentary) are individually packaged, labeled and provided **NON-STERILE**.
- Each bone screw should be measured for its length (using the measurement scale on the screw caddy) and placed in its designated hole in the screw caddy.
- Before each use (surgical procedure), the screw caddy needs to be completely assembled. The bone screws (in the screw caddy) need to be **CLEANED** and **STERILIZED** along with the surgical instruments in the instrument tray in accordance with the recommended procedures (see below).

**CLEANING**

Remove and discard any plastic caps or tip protectors before cleaning and sterilizing implants and instruments. Cleaning by hand rather than by mechanical cleaning will prolong the life of the implants and instruments. Clean all crevices, flutes, and cannulations of all debris, using a soft bristle brush or cleaning stylet. Remove all traces of blood or other residues immediately. Do not allow these to dry. Implants and instruments should be cleaned while submerged in warm water with an appropriate neutral pH detergent. Always follow the manufacturer’s instructions when preparing and using detergents. Do not use steel brushes as they can accelerate wear and corrosion of the implants and instruments. Rinse implants and instruments thoroughly with distilled water. Dry implants and instruments immediately after cleaning.

**STERILIZATION**

Recommended parameters for steam sterilization are as follows:

Cycle	Temperature	Minimum Exposure Time
Vacuum	270° F/ 132° C	4 minutes
Gravity	250° F/ 121° C	30 minutes

Recommended dry time is 30 minutes.

Parameters may vary based on manufacturer, installation or maintenance of sterilization equipment. On-going testing must be performed by the user to confirm inactivation of all forms of microorganisms.

Sterilizing in liquid solutions is not recommended. Do not sterilize at temperatures greater than 275° F/ 135° C.

**CAUTION**

Federal Law (USA) restricts this device to sale by or on the order of a physician.

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All rights reserved. Printed in the U.S.A
ArthroSurface, Inc.
MT Promedt Consulting GmbH
ArthroSurface, Inc. U.S.© 2019 ArthroSurface, Inc.