**Unicompartmental Knee  
Resurfacing Prosthesis (UniCAP™)  
Instructions for Use**

## Unicompartmental Knee Resurfacing Prosthesis (UniCAP™) Instructions for Use

### Description

The Unicompartmental Knee Resurfacing Prosthesis (UniCAP™) incorporates a low-profile femoral articular component that mates to a taper post via a taper interlock. The femoral component articulates against an all-polyethylene tibial component. The UniCAP™ implants allow resurfacing of the compartment utilizing the undisturbed compartmental structures and soft-tissues.

### Materials

Femoral Components	
Articular Resurfacing Component:	Cobalt-Chromium-Molybdenum alloy (Co-Cr-Mo)
Surface Coating:	Titanium (CP Ti)
Taper Post:	Titanium Alloy (Ti-6Al-4V)
Tibial Components	
Ultra-High Molecular Weight Polyethylene (UHMWPE)	

### Indications for Use

Partial replacement of the articulating surfaces of the knee when only one side of the joint is affected due to the compartmental primary degenerative or post-traumatic degenerative disease, previous tibial condyle or plateau fractures, deformity or revision of previous arthroplasty. This device is intended to be used with bone cement.

#### Patient selection factors to be considered include:

1. Patient's need to obtain pain relief and improve function is significant.
2. Patient's joint stability is good, with intact ACL, and limited mechanical axis deformity.
3. Patient's overall well-being is good, including an ability and willingness to follow instructions and comply with activity restrictions.

### Contraindications

#### Absolute contraindications include:

1. Infection, sepsis, or osteomyelitis.

2. Patients that have a known sensitivity to materials typically used in orthopedic prosthetic devices or bone cements.

**Relative contraindications include:**

1. Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
2. Osteoporosis.
3. Metabolic disorders which may impair the formation or healing of bone.
4. Infections at remote sites which may spread to the implant site.
5. Rapid joint destruction or bone resorption visible on roentgenogram.
6. Chronic instability or deficient soft tissues and other support structures.
7. Vascular or muscular insufficiency.
8. Inadequate skin, musculotendinous or neurovascular system status.

**Warnings**

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

The all-polyethylene tibial component is intended to be used in conjunction with activity restrictions.

When taking readings of articular surfaces, care should be taken to ensure that the distal tip of contact probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on contact probe to slightly indent articular surface at chosen points, this will ensure selected implant will be recessed just below articular surface at margins of implant.

When placing implant, carefully trim articular cartilage debris around margin of implant. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the taper post and implant, carefully clean taper post with provided instruments. All drilling or reaming should be done with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Ensure that care is taken to obtain complete and uniform bone cement coverage at implant site. Unsupported components or unevenly supported components may result in implant failure.

Accepted practices in post-operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post-operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

#### **Precautions**

The Unicompartamental Knee Resurfacing Prosthesis (UniCAP™) implant is intended to be fitted and installed with the Unicompartamental Knee Resurfacing Prosthesis (UniCAP™) instrument set. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The instrument set should be regularly inspected for any signs of wear or damage. Do not reuse implants. Reuse of single use devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device.

#### **Possible Adverse Effects**

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or post-operative bone fracture.
9. Post-operative pain or incomplete resolution of preoperative symptoms.

10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient peroneal palsy.

#### **Sterility**

Metallic prosthetic components are sterilized by exposure to gamma irradiation. Non-metallic prosthetic components are sterilized by gas plasma sterilization. Do not resterilize any components. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date.

**Caution: United States Federal law restricts this device to sale by or on the order of a physician.**

Manufactured for ArthroSurface<sup>®</sup>, Inc. by:

STD Med, Inc  
75 Mill Street  
Stoughton, MA 02072

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1. The tibial component may be placed under arthroscopic visualization to ease preparation of the tibial defect. With the knee at 90 degrees of flexion, and working through an anteromedial portal, place a **Tibial Template** central to the damaged area of the tibial plateau surface. The underside of the **Tibial Template** should be in contact with the surface of the tibia. Select the **Tibial Template** that best matches the A/P and M/L surface curvatures of the tibia. Attach the **Tibial Template** to the **Tibial Drill Guide** and **Bullet**. Place the guide so the arm of the **Tibial Template** is parallel to the tibial plateau.

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
2. Drill the **2.0mm Drill Tipped Pin** through the central axis of the **Tibial Drill Guide** until it reaches the center of the **Tibial Template**. Use a small closed curette to “catch” the tip of the **Pin** to prevent drilling into the femur. Care must be taken so that excessive torque is not applied to the **Drill Guide** which may cause the pin to miss the target. Confirm that there is a minimum of 5mm of bone from the edge of the **Template** to front of the tibia to avoid breaking through the anterior tibia during reaming. Bring the knee into extension and probe the underside of the meniscus to ensure complete visualization and proper placement of the **Template**. Remove the **Tibial Drill Guide** and **Bullet**.

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
3. Drive the **Tibial Pilot Drill** over the **2.0mm Drill Tipped Pin** until it reaches the center of the templated area. Drilling should stop before the larger diameter tip of the **Pilot Drill** breaches the tibial plateau. Remove the **Tibial Pilot Drill** and **2.0mm Pin**.

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4. Advance the **Introducer** into the prepared tibial tunnel. The proximal tip of the **Introducer** should be flush with the tibial plateau. Begin to advance the threaded **Blade Stop** over the **Introducer** until it begins to screw into the bone. Remove **Introducer** and continue to advance **Blade Stop** until it is 2/3rds into the tunnel.
  - a. Remove **Driver** handle and reinsert **Introducer** and **Driver** and continue to advance as one unit. Stop when the tip of the **Introducer** is flush with the tibial plateau. Confirm that the laser mark on the **Introducer** is in-line with the laser mark in the slotted window of the **Blade Stop Driver**.
  - b. **Blade Stop** is at the correct depth when the tip of the **Introducer** is flush with the tibial plateau and the laser mark lines on the **Driver** and **Introducer** are aligned in the slotted window. Remove the **Blade Stop Driver** and **Introducer**.

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
5. Place the **Cutting Blade** into the **Blade Holder** with the long slot facing posteriorly to the joint. Introduce the **Cutting Blade** into the portal. Advance the **Blade Drive Shaft** into the tibial tunnel until it is visible in the joint. Push the tip of the **Drive Shaft** through the center of the **Cutting Blade**.

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6. To lock the Tibial Cutting system, push the sheath in an upward motion and rotate 90 degrees so the Lock Indicator on **Blade Drive Shaft** is positioned over the Dowel Pin and Laser Mark line. Release to lock **Cutting Blade** into position. Attach the

Powered Drill to the laser mark indicated on the distal end of the **Drive Shaft**.

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7. Using a Powered Drill, begin rotation counterclockwise to normalize the blade to the tibial plateau. This will help ensure even cutting engagement of the **Cutting Blade** into the plateau. Care should be taken to avoid the meniscus. Begin to prepare the inlay implant socket using a clockwise blade rotation. Drilling is complete when the **Cutting Blade** comes in contact and stops on the proximal end of the **Blade Stop**.

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
8. With **Cutting Blade** remaining in place, insert the appropriate sized **Sizing Trial** based on the mapping determined with the

**Tibial Template**. Confirm fit at anterior/posterior and medial/lateral margins of the **Sizing Trial**. If the **Trial** is proud at the margins, use the **Blade Stop Wrench** to rotate the **Blade Stop** clockwise. Each 90 degree turn of the **Blade Stop** with **Blade Stop Wrench** lowers the **Blade Stop** 1mm. Reattach the Powered Drill and re-ream to the new depth. This will lower the implant depth in the tibial socket.

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9. With the **Sizing Trial** set at the appropriate height, begin removing the instrumentation. Raise the **Drive Shaft** so the **Cutting Blade** can be grasped for removal. To unlock and remove the **Cutting Blade** push the sheath in an upward motion and rotate 90 degrees counter clockwise so the Unlock Indicator

on **Blade Drive Shaft** is positioned over the Dowel Pin. Release to unlock **Cutting Blade**. Push **Drive Shaft** upwards to free it from the **Cutting Blade** and turn the **Drive Shaft** 90 degrees. Pull distally on the **Drive Shaft** and remove from the tibial tunnel. Use a grasper to remove the **Cutting Blade**.

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
10. Open **Tibial Component Kit**. Using the **Suture Retriever** capture the suture and pull suture through and out of the distal tibial drill hole. Introduce the **Tibial Component** into the tibial socket using the **Delivery Tool**.

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
11. Advance the **Slot Driver** into the tibial tunnel to rotate the **Tibial Component** (via the distal slot on the bottom of the component) to its optimal orientation if needed.

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
12. Implant **Tibial Component** with cement using the **Arthrosurface® Cement Ejector**. When discharging cement, allow the back pressure from the cement extraction to lift the

implant up 2mm and then continue to back fill the tibial tunnel with cement. Utilize the **Tibial Template** through the portal to apply downward pressure onto the **Tibial Component** to seat it in its final position. This will allow for optimal cement integration.

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
#### **Cement Ejector Assembly**

- Mix low viscosity cement according to manufacturers' directions.
- Place cement into the **Delivery Syringe**. Remove **Funnel** when **Delivery Syringe** is full.
- Insert **Plunger** into **Delivery Syringe**. Insert into **Ejector Handle**.
- Attach **Drive Rod** onto **Ejector Handle**.
- Place **Threaded Sheath** into tibial tunnel to prevent cement extrusion.
- A powered drill with Jacobs chuck is used to advance the **Drive Rod**.
- Deliver cement to undersurface of implant through tibial tunnel and draw **Cement Ejector** retrograde as tibial tunnel is filled.

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
1. With knee at 90 degrees of flexion and working through an anteromedial incision, determine the anterior/posterior and medial/lateral curvatures of the condyle using the **Templates**.

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
2. Select the **Sizing Jig** based on the determined anterior/posterior curvature. Position and secure the **Sizing Jig** on the condyle using the **2.0mm Short Pins**.

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

3. Center the **Bushing** into the **Sizing Jig**, beginning with the most inferior position. Drill the **Short Threaded Pin** through the **Bushing** lumen into bone, stopping when the laser mark line on the **2.0mm Threaded Pin** is level with the proximal surface of the bushing. Remove the **Bushing**.

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
4. Advance the **17.5mm Reamer** over the **Threaded Pin** until it contacts the stop in the slotted proximal window of the **Reamer**. Repeat steps 3 and 4 for the remaining two positions in the **Sizing Jig**.

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
5. Select the appropriately sized **Femoral Trial**, attach to **Sizing Trail Handle** and place into position. Confirm fit at anterior/posterior margins and medial/lateral margins of the **Femoral Trial**. **Femoral Trial** should be seated so the edges are slightly recessed relative to the surrounding cartilage.

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
6. Confirm position of **Femoral Trial** before preparing the pilot hole for the **Taper Post**.

Advance the **4.7mm Pilot Drill** through the **Sizing Trial Handle** until the laser mark is flush with the end of the handle and leave it in position. Replace the **2.0mm Pin** through the **4.7mm Pilot Drill**.

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
7. Remove **Sizing Trial Handle** and **4.7mm Pilot Drill**, leaving **2.0mm Pin** in position. Place **Step Drill** over the **2.0mm Pin** and reposition **Sizing Trial Handle**. Advance **Step Drill** until the laser mark is flush with the end of the handle.

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
8. Remove Sizing Trial Handle and Step Drill, leaving 2.0mm Pin in position. Feed proximal end of Tap into distal opening of Sizing Trial Handle, place unit over the 2.0mm Pin and reposition Sizing Trial Handle. Couple the Tap Handle to the Tap and advance the Tap until the laser mark is flush with the end of the handle.

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9. Insert the **Taper Post** into the **Sizing Trial Handle**. Attach the **Handle** and **Taper Post** assembly onto the **Femoral Sizing Trial**. Insert the **Hex Driver** into the **Handle** and advance the **Taper Post** into the bone. Stop advancing the **Hex Driver** when the raised stop on the **Driver Shaft** contacts the top of the **Sizing Trial Handle** and the **Femoral Sizing Trial** is flush with the surrounding cartilage. Remove the **Sizing Trial**.

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10. Use the **Final Placement Gauge** in combination with the **Femoral Trial** to confirm the proper depth placement of the **Taper Post**. With the **Femoral Trial** in position, and with the **Final Placement Gauge** coupled to the **Taper Post**. A very small separation (less than .5mm) should be seen between the components. If no gap is visible, the **Taper Post** is set too deep and should be raised. If a large gap is visible, the **Taper Post** is set too shallow and should be lowered.

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
11. Apply pea-sized balls of bone cement to underside of **Femoral Component**. Position **Femoral Component**. Use a slight tap on the **Impactor** to mate **Femoral Component** to **Taper Post**.

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**R<sub>x</sub>ONLY**

PN K001-2002 REV D

This product is covered by one or more of U.S. Patent Nos. 6,520,964;  
6,610,067; 6,679,917; other patents and other patents pending. HemiCAP®  
is a trademark of ArthroSurface, Inc. U.S.

© 2018 ArthroSurface, Inc. All rights reserved. Printed in the U.S.A.

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RSQR Ltd.