WristMotion® Hemiarthroplasty System  
Instructions for Use

---

## **WristMotion® Hemiarthroplasty System Instructions for Use**

---

### **Description**

The Arthrosurface WristMotion® Hemiarthroplasty System consists of a contoured capitate articular implant designed to articulate with the natural radius bone, a taper post and set of instruments used for implant site preparation and delivery. The capitate articular components are manufactured using implant grade cobalt chrome alloy and is offered in two diameters, and four articular radii. The taper post is manufactured using implant grade titanium alloy and is offered in one fixed size designed to work with all capitate articular implants.

### **Materials**

Articular Component:

Cobalt-Chrome Alloy (Co-Cr-Mo)

Surface Coating: Titanium (CP Ti)

Taper Post Component: Titanium Alloy (Ti-6Al-4V)

### **Indications for Use**

The Arthrosurface WristMotion® Hemiarthroplasty System is indicated for use as a partial replacement of wrist joint(s) disabled by pain, deformity and/or limited motion caused by:

1. Non-inflammatory degenerative joint disease including osteoarthritis, traumatic arthritis and avascular necrosis.
2. Rheumatoid arthritis.
3. Revision where other devices or treatments have failed.
4. Scapholunate Advanced Collapse (SLAC) and other functional deformities.
5. Trauma, including fractures of the carpal bones.

The device is a single-use implant intended to be used with bone cement.

**Patient Selection Factors to be Considered Include:**

1. Need to obtain pain relief and improve function.
2. Patient age as a potential for early-age-revision of total joint arthroplasty or arthrodesis.
3. Patient overall well-being, including ability and willingness to follow pre and post-operative treatment regimen.
4. Failure of previous conservative treatment options in correcting deformity and achieving pain relief.

**Contraindications**

**Absolute contraindications include:**

1. Significant bone demineralization or inadequate bone stock.
2. Inadequate skin, musculotendinous or neurovascular system status.
3. Infection, sepsis and osteomyelitis.
4. Patients that have a known sensitivity to cobalt-chrome and/or titanium alloys typically used in prosthetic devices.

**Relative contraindications include:**

1. Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
2. Osteoporosis.
3. Metabolic disorders which may impair the formation or healing of bone.
4. Infections at remote sites which may spread to the implant site.
5. Rapid joint destruction or bone resorption visible on roentgenogram.
6. Chronic instability or deficient soft tissues and other support structures.
7. Vascular or muscular insufficiency.
8. Absent or insufficient wrist extensor tendons.

### **Warnings and Precautions**

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When placing implant, carefully trim articular cartilage debris around margin of implant. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the taper post and articular component, carefully clean taper with provided instruments. All drilling or reaming should be done at slowest speeds possible with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Accepted practices in postoperative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to postoperative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

All surgical implants are subjected to repeated stresses that can result in failure. The use of an implant should be avoided if excessive loading cannot be prevented at or near the implant site. No other metallic or non-metallic implantable devices are to be used in conjunction with the Arthrosurface WristMotion® Hemiarthroplasty System at the implant site. Doing so may compromise implant performance and patient safety.

These implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. Their safety in the MR

environment is unknown. Scanning a patient who has this device may result in patient injury.

#### **Precautions**

These implants are intended to be fitted and installed with the associated instruments. Use of instruments from other systems may result in improper implant selection, fitting and placement, which could result in implant failure or poor clinical outcome. The instruments should be regularly inspected for any signs of wear or damage.

Surgeon or Physician should discuss general risks and potential complications associated with this and any surgical procedure with the patient prior to patient consent.

#### **Possible Adverse Effects and Complications**

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.

8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.
13. Embolism.

**Sterility**

Implants and single-use disposable instruments are provided STERILE. Metallic implant components are sterilized by exposure to gamma radiation. Do not resterilize. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date. Do not reuse implants or single-use disposable instruments. Reuse of these devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device(s).

**Caution**

United States Federal Law restricts this device to sale by or on the order of a physician.

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
- Using the **Coronal Radius Template** and the **Sagittal Radius Template** determine the appropriate curvatures of the articular surfaces of the distal radius. This combination corresponds with a prescribed capitate articular implant curvature. Use the decision matrix below to

determine the indicated curvature. Note: the **12 mm Articular Component** is only available with a sagittal curvature of 17 mm.

Coronal		Sagittal	
Radius Curvature (Measured)	Capitate Curvature (Indicated)	Radius Curvature (Measured)	Capitate Curvature (Indicated)
If <b>34 mm</b> Measured	Use <b>22 mm</b>	If <b>18 mm</b> Measured	Use <b>17 mm</b>
If <b>56 mm</b> Measured	Use <b>35 mm</b>	If <b>25 mm</b> Measured	Use <b>23 mm</b>

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7. With the **Drill Guide**, locate the capitate's axis normal to the articular surface. Use the **Drill Guide** to determine which capitate articular component diameter matches the surface diameter. Place **1.5 mm Guide Pin** into a cannulated power drill. Advance the **1.5 mm Guide Pin** into the bone making sure that is central to the articular surface. It is important to verify that the **Drill Guide** is seated on the surface such that all 4 points of contact are established on the articular surface. A normal axis and correct diameter measurement are necessary for proper implant fit.

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11. Choose the appropriate **Capitate Reamer** based on the diameter measured by the **Drill Guide** in Step 7. Drive **Capitate Reamer** over **1.5 mm Guide**

**Pin** until it contacts the top surface on the **Taper Post** (Use lavage during drilling to prevent possible tissue damage from heating effects). Make sure not to bend the **1.5 mm Guide Pin** during drilling as it may result in **Articular Component** misalignment.

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HemiCAP is a trademark of ArthroSurface, Inc.
© 2020 ArthroSurface, Inc. All rights reserved.
ArthroSurface, Inc.
MT Promedt Consulting GmbH
Printed in the U.S.A.