## **Contoured Articular Prosthetic (CAP) Resurfacing System Instruments**

### **Description**

The Contoured Articular Prosthetic HemiCAP<sup>®</sup> instruments are designed to be used in the sizing, location, and delivery of HemiCAP<sup>®</sup> fixation components and articular components. These instruments are designed for repeated use, with proper care and handling.

### **Instructions for Use**

Instructions for use for each of the HemiCAP<sup>®</sup> instruments is provided in the HemiCAP<sup>®</sup> implant *Instructions for Use* package insert.

Check Instrument Set and implant labeling to ensure product match prior to usage.

### **Warnings and Precautions**

This product is provided **NON-STERILE**. The product must be properly cleaned and sterilized before each use. Remove and discard any plastic caps or tip protectors before cleaning and sterilizing the instruments. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing the HemiCAP<sup>®</sup> procedure. The HemiCAP<sup>®</sup> instrument set should be regularly inspected for any signs of wear or damage.

### **Cleaning**

Cleaning instruments by hand rather than by mechanical cleaning will prolong the life of the instruments. Clean all crevices, flutes, and cannulations of all debris, using a soft bristle brush or cleaning stylet. Remove all traces of blood or other residues immediately. Do not allow these to dry. Instruments should be cleaned for one minute while submerged in a warm water bath filled with an appropriate neutral detergent. All detergents should be prepared per use-dilution and temperature recommendations of the manufacturer. Do not use

steel brushes as they can accelerate wear and corrosion of the instruments. Rinse instruments thoroughly for one minute with distilled water. Lumens and other hard to reach locations should be thoroughly flushed. Dry instruments immediately after cleaning.

### **Sterilization**

Recommended parameters for steam sterilization are as follows:

<u>Cycle</u>	<u>Temperature</u>	<u>Minimum Exposure Time</u>
Vacuum	270° F/ 132° C	20 minutes
Gravity	250° F/ 121° C	60 minutes

Recommended dry time is 30 minutes.

Parameters may vary based on manufacturer, installation or maintenance of sterilization equipment. On-going testing must be performed by the user to confirm inactivation of all forms of microorganisms.

Instrument Trays should be processed in double wrapped configuration using an FDA (or applicable regulatory agency) cleared sterilization wrap.

Sterilizing in liquid solutions is not recommended. Do not sterilize at temperatures greater than 275° F/ 135° C.

### **Inspection**

All instruments should be routinely evaluated for visible wear and damage. Instruments which function with the aid of a central working axis should be evaluated by performing a functional test to ensure no occlusions exist within the instrumental lumens. Each of these instruments should be tested by sliding the instrument over the designated central working axis piece per the *Instructions for Use* package insert. Should the instruments fail to engage properly, they should be deemed unfit for use.

**Note:** ArthroSurface does not define the maximum number of uses appropriate for reusable instruments. The useful life of these devices depends on many factors, including the method and duration of each use and handling between uses.

**CAUTION**

United States Federal Law restricts this device to sale by or on the order of a physician.

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© 2020 ArthroSurface, Inc. All rights reserved. Printed in U.S.A.
Arthrosurface, Inc. 28 Forge Parkway, Franklin, MA 02038,
ArthroSurface, Inc.
MT Promedt Consulting GmbH