---

**Description**

The Patello-Femoral Wave Arthroplasty Systems incorporate a distal femoral trochlear surface articular component that mates to a taper post via a taper interlock, and an all-polyethylene patella component. The prosthesis is intended to be used in cemented arthroplasty.

**Materials**

Femoral Component:

Cobalt-Chromium Alloy (Co-Cr-Mo)

Undersurface Coating: Titanium (CP Ti)

Taper Post: Titanium Alloy (Ti-6Al-4V)

Patella Component:

Ultra-High-Molecular Weight Polyethylene (UHMWPE)

**Indications**

Intended to be used in cemented arthroplasty in patients with osteoarthritis limited to the distal patello-femoral joint, patients with a history of patellar dislocation or patellar fracture, and those patients with failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release, etc.) where pain, deformity or dysfunction persists.

**Patient selection factors to be considered include:**

- 1) Patient's need to obtain pain relief and improve function is significant;
- 2) Patient's tibio-femoral joint is substantially normal;
- 3) Patient exhibits no significant mechanical axis deformity;
- 4) Patient's menisci and cruciates are intact with good joint stability, and good range of motion; and
- 5) Patient's overall well-being is good, including the ability and willingness to follow instructions and comply with activity restrictions.

**Contraindications****Absolute contraindications include:**

- 1) Defects that are not localized.
- 2) Inflammatory degenerative joint disease, rheumatoid arthritis, infection, sepsis, or osteomyelitis.
- 3) Patients that have a known sensitivity to materials typically used in orthopedic prosthetic devices or bone cements.

**Relative contraindications include:**

- 1) Uncooperative patient or patient incapable of following pre-operative and post-operative instructions.
- 2) Metabolic disorders, which may impair the formation or healing of bone; osteoporosis.

- 3) Infections at remote sites, which may spread to the implant site.
- 4) Rapid joint destruction or bone resorption visible on roentgenogram.
- 5) Chronic instability or deficient soft tissues and other support structures.
- 6) Vascular or muscular insufficiency.
- 7) Inadequate skin, musculotendinous or neurovascular system status.

#### **Warnings**

Improper selection, placement, positioning, alignment, and fixation of the implant components may reduce the service life of the prosthetic components. Inadequate preparation and cleaning of the implant components mating surfaces may result in improper fixation of the device. Improper handling of the implants can produce scratches, nicks or dents that may have adverse clinical effects on mating joint surfaces. Do not modify implants. The surgeon shall be thoroughly familiar with the implants, instruments, and surgical technique prior to performing surgery.

When defining offsets of articular surfaces, care should be taken to ensure that instruments are properly aligned and mated with taper in Taper Post. Visually confirm distal tip of Contact Probe is making contact on articular surfaces and free from any soft tissue structures to ensure accuracy. Use light pressure on Contact Probe to slightly indent articular surface at each mapping point, ensuring that the selected implant will be flush or slightly recessed with the articular surface.

Prior to placing implant, carefully trim articular cartilage debris around prepared margin. Remove bone particles and lavage thoroughly. To ensure mechanical interlock of the Taper Post and Femoral Component, carefully clean Taper Post taper with provided instruments. All drilling or reaming should be done with vigorous lavage to minimize heat effects to adjacent bone and cartilage tissues.

Ensure that care is taken to obtain complete and uniform bone cement coverage at implant site. Unsupported components or unevenly supported components may result in implant failure.

Accepted practices in post operative care should be used. The patient is to be instructed and monitored to ensure a reasonable degree of compliance to post operative instructions and activity restrictions. Excessive activity, impact, and weight gain have been implicated in the reduction of the benefit and service life of prosthetic devices.

The Patello-Femoral implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of The Patello-Femoral knee implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### **Precautions**

Patello-Femoral implants are intended to be fitted and installed with the appropriate Patello-Femoral instrument sets. Use of instruments from other systems may result in improper implant selection, fitting, and placement which could result in implant failure or poor clinical outcome. The instrument set should be regularly inspected for any signs of wear or damage. Do not reuse implants or disposable instruments. Reuse of single use devices can increase the risk of patient infection and can compromise service life and other performance attributes of the device.

#### **Possible Adverse Effects**

1. Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions. Particulate wear debris and mild tissue discoloration from metallic components have been noted in other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
2. Infection or allergic reaction.
3. Loosening, migration or loss of fixation of implant.
4. Fretting and crevice corrosion can occur at the interface between the implant components.
5. Fatigue fracture of the implants as a result of bone resorption around the implant components.
6. Wear and damage to the implant articulating surface.
7. Wear and damage to the adjacent and opposed articular cartilage surfaces or soft tissue support structures.
8. Intraoperative or postoperative bone fracture.
9. Postoperative pain or incomplete resolution of preoperative symptoms.
10. Periarticular calcification or ossification, with or without impediment of joint mobility.
11. Incomplete range of motion due to improper selection or positioning of components.
12. Transient nerve palsy.

### **Sterility**

Metallic prosthetic components are sterilized by exposure to gamma irradiation. Non-metallic prosthetic components are sterilized by gas plasma sterilization. Do not resterilize any components. Do not use components if packaging is opened or damaged. Do not use components if beyond expiration date.

Caution: United States Federal law restricts this device to sale by or on the order of a physician.

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7. Assemble the **Outer Reamer** into the **Guide Bushing**. Secure **Guide Bushing** into the **Guide Block**. Advance **Outer Reamer** into the superior **Guide Block** opening and into bone until depth mark is reached. Repeat for the inferior **Guide Block** opening. Remove **Guide Block**. If using the **Kahuna Femoral Component**, repeat the superior ream only, using the disposable **Kahuna Reamer**.

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14. Use the **Final Placement Gauge** in combination with the **Sizing Trial** to confirm the proper depth placement of the **Taper Post**. With the **Sizing Trial** in position, and with the **Final Placement Gauge** coupled to the **Taper Post**. A very small separation (less than .5mm) should be seen between the components. If no gap is visible, the **Taper Post** is set too deep and should be raised. If a large gap is visible, the **Taper Post** is set too shallow and should be lowered. Be sure not to impact the **Final Placement Gauge** onto the **Taper Post**.

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
15. Prior to placing the **Femoral Component** on the **Implant Holder (suction)** make sure that sufficient suction is present to hold the device on the distal suction cup. Align the **Femoral Component** on the **Implant Holder** with the medial etch mark facing the medial aspect of the knee and lateral mark facing the lateral plane. Apply a small amount of bone cement to the backside of the **Femoral Component**. Insert **Femoral Component** into taper of the Taper Post.

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
17. Alternative fixation is to pre-assemble the **Threadless Stud** to the **Femoral Component**. Be sure to protect the articular face of the **Femoral Component**, use slight impaction with the mallet to seat the Morse taper of the **Threadless Stud** onto the **Femoral Component**. Apply a small amount of bone cement to the backside of the **Femoral Component**. Using the **Impactor**, firmly mallet the **Femoral Component** until completely seated.

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
1. With knee in extension, locate the **Alignment Guide** so that the pin fits into the **Taper Post**. While observing range of motion, identify target placement of the **Patella Component** using the pointer on the **Alignment Guide** to transfer **Taper Post** central axis. Use slight pressure against the patella so that the pointer on the **Alignment Guide** creates an indentation on distal patella surface.

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
2. Place the **Drill Guide** so that its central axis passes through the indentation on the patella surface. Drill the **Guide Pin** through until it engages the opposite cortex of the patella. *(It is important to verify that the **Drill Guide** is seated on the curved surface such that all 4 points of contact are established on the articular surface. A normal axis is necessary for proper implant fit).*

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
11. Using the **Patella Clamp**, place the *Anatomic or Button* contacting surface against the **Patella Component** and the opposed surface in-line behind the patella. Tighten the **Patella Clamp** until the **Patella Component** is firmly seated in the prepared socket. Leave the **Patella Clamp** in place while the bone cement adequately cures. Remove the **Patella Clamp** and clean out any remaining exposed cement.

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CE  
0459

R<sub>x</sub> ONLY

This product is covered by one or more of U.S. Patent  
Nos. 6,520,964; 6,610,067; 6,679,917; other patents and  
other patents pending. HemiCAP® is a trademark of  
Arthrosurface, Inc. U.S.

© 2022 Arthrosurface, Inc. All rights reserved.  
Printed in the U.S.A.

PN 2001-2110 REV I

Patello-Femoral Wave Arthroplasty Systems  
Instructions for Use

*PF Wave® and WaveKahuna® are marketed and sold under the trademarks PFXL™ and PF Kahuna™ in the EU and UK.*

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MT Promedt Consulting GmbH